

Tracoe Softpad – REF 969

Dressings

Product Description

There are many essential products that play an important role in tracheostomy care. Our assortment of essentials includes neckbands, dressings, products for care & cleaning as well as clothing & protectors and accessories. They fulfill a wide range of user needs and preferences, are easy to use, comfortable, discreet, offer security, and are functional.



General Information		Material Specifications	
Product line	Dressings	Material Compress	Polyurethane foam
Legal manufacturer	Tracoe Medical GmbH Reichelsheimer Str. 1/3 55268 Nieder-Olm Germany	Biocompatibility	ISO 10992-5 & ISO 10993-10
Country of origin	Germany		
Handling		Storage	
Shelf life	3 years		
Maximum Useful Life	Change after cannula cleaning or because of heavy secretions	Certification	
Packaging		Product classification acc. MDR	I
Packaging unit	20 pcs	CE Marking	 0 1 2 3
Dimensions (W x H x L)	145 x 100 x 105 mm	GMDN	46854
Material sterile packaging	Paper poly pouch	Quality system certified to	EN ISO 13485
Material outer packaging:	Cardboard box	FDA registration	Exempt/E652368

Product Specifications				
Product Name	REF	Dimensions (W x L x H)	Dimensions of Hole	Slit
Tracoe Softpad	REF 969	3.6 x 5.8 x 0.5 cm	10 mm	4 small slits around Hole



Tracoe Medical GmbH,
Reichelsheimer Str. 1 / 3,
55268 Nieder-Olm
Germany



Breathing-Speaking-Living

Atos Medical AB
Hyllie Boulevard 17
215 32 Malmö
Sweden
www.atosmedical.com