

Tracoe Twist – 501-XXX

Tracoe Twist spare inner cannulas

Product Description

The Tracoe Twist tracheostomy tubes are used as ventilation and therapy tubes with inner cannulas in hospitals, rehabilitation, and community settings. A special feature is the curved neck flange that moves in both the vertical and horizontal axes, following the movements of the patient's head. The thin tube walls of inner and outer cannulas allow for a good inside - outside diameter ratio. Tracoe Twist and Twist Plus tracheostomy tubes are made of polyurethane (PU). The low-pressure cuff is made of medical grade PVC.



Scope of Supply

- 30 inner cannulas with 15 mm connector

General Information	
Product line	Twist
Legal manufacturer	Tracoe Medical GmbH Reichelsheimer Str. 1/3 55268 Nieder-Olm, Germany
Country of origin	Germany
Colour Coding	n. A.
Handling	
Shelf life	5 years, from date of production
Maximum Useful Life	29 days, from first use

Packaging	
Packaging unit	30 units
Weight	See Table
Dimensions (W x H x L)	320 x 440 x 0,07
Material outer packaging:	Cardboard box
Material sterile packaging:	Sterile bag with medical paper and transparent heatsealable PET/PP foil

Material Specifications	
Material Outer cannula	Polyurethane
Material Neck Flange	n. A.
Material Cuff	n. A.
Biocompatibility	EN ISO 10993-1

Storage	

Certification	
Product classification acc. MDR	IIb
CE Marking	
GMDN	35404
Quality system certified to	EN ISO 13485
FDA registration	Clearance Number K961449

Product Range

Size	Order Number	ID-IC	BA	Gross Weight	Global Trade Item Number
		mm	∅	g	
06	REF 501-06-XXX	6,3	95	209	04035324027965
07	REF 501-07-XXX	7,3	95	212	04035324027958
08	REF 501-08-XXX	8,0	90	211	04035324027866
09	REF 501-09-XXX	9,0	90	209	04035324027972

Legend

ID-IC Inside diameter at bottom of inner cannula

BA Bending angle



Tracoe Medical GmbH
Reichelsheimer Str. 1 / 3
55268 Nieder-Olm
Germany



Breathing-Speaking-Living

Atos Medical AB
Hyllie Boulevard 17
215 32 Malmö
Sweden
www.atosmedical.com